



Inbios International, Inc.  
562 1<sup>st</sup> Ave. South  
Suite 600  
Seattle, WA 98104  
206-344-5821 (phone)  
206-344-5823 (fax)  
[www.inbios.com](http://www.inbios.com)

## Product Specification Sheet

Product	Monoclonal antibody specific to <i>Plasmodium falciparum</i> (and other <i>Plasmodium</i> species)
Catalog Number	INB1611
Description	Antibody specific to <i>Plasmodium falciparum</i> , <i>vivax</i> , and <i>ovale</i> LDH
Specificity	LDH
Clone	1017.161
Host Animal	Mouse
Format	Purified, liquid
Isotype	IgG1
Buffer	PBS pH 7.4
Application	Suitable for use in ELISA and Rapid Immunoassays. Each laboratory should determine an optimum working titer for their own application. Other formats have not been tested.
Storage	Store at -20°C. Avoid multiple freeze/thaw cycles. Prepare working dilution only prior to immediate use.

**For Research Use Only. Not for use in Diagnostic Procedures.**

**Important Note: Centrifuge before opening to ensure complete recovery of vial contents.**