



Inbios International, Inc.
562 1st Ave. South
Suite 600
Seattle, WA 98104
206-344-5821 (phone)
206-344-5823 (fax)
www.inbios.com

Product Specification Sheet

| | |
|----------------|--|
| Product | Monoclonal antibody specific to <i>Plasmodium falciparum</i> |
| Catalog Number | INB1760 |
| Description | Antibody specific to <i>Plasmodium falciparum</i> LDH |
| Specificity | LDH |
| Clone | 915.7601 |
| Host Animal | Mouse |
| Format | Purified, liquid |
| Isotype | IgG1 |
| Buffer | PBS pH 7.4 |
| Application | Suitable for use in ELISA and Rapid Immunoassays. Each laboratory should determine an optimum working titer for their own application. Other formats have not been tested. |
| Storage | Store at -20°C. Avoid multiple freeze/thaw cycles. Prepare working dilution only prior to immediate use. |

For Research Use Only. Not for use in Diagnostic Procedures.

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.