



Inbios International, Inc.  
562 1<sup>st</sup> Ave. South  
Suite 600  
Seattle, WA 98104  
206-344-5821 (phone)  
206-344-5823 (fax)  
[www.inbios.com](http://www.inbios.com)

## Product Specification Sheet

Product	Monoclonal antibody specific to <i>Shigella sonnei</i>
Catalog Number	INB3100
Description	Antibody specific to <i>Shigella sonnei</i> lipopolysaccharide
Specificity	LPS
Clone	610.9100
Host Animal	Mouse
Format	Purified, liquid
Isotype	IgG1
Buffer	PBS pH 7.4
Application	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for their own application. Other formats have not been tested.
Storage	Store at -20°C. Avoid multiple freeze/thaw cycles. Prepare working dilution only prior to immediate use.

**For Research Use Only. Not for use in Diagnostic Procedures.**

**Important Note: Centrifuge before opening to ensure complete recovery of vial contents.**