



Inbios International, Inc.
562 1st Ave. South
Suite 600
Seattle, WA 98104
206-344-5821 (phone)
206-344-5823 (fax)
www.inbios.com

Product Specification Sheet

Product	Monoclonal antibody specific to Dengue Virus
Catalog Number	INB5283
Description	Antibody specific to Dengue Virus nonstructural protein 1 (NS1)
Specificity	DENV1,2, 3, 4
Clone	528.1133
Host Animal	Mouse
Format	Purified, liquid
Isotype	IgG1
Buffer	10mM TrisCl 50mM NaCl pH 8.0
Application	Suitable for use in ELISA and Rapid Immunoassays. Each laboratory should determine an optimum working titer for their own application. Other formats have not been tested.
Storage	Store at -20°C. Avoid multiple freeze/thaw cycles. Prepare working dilution only prior to immediate use.

For Research Use Only. Not for use in Diagnostic Procedures.

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.