



InBios International, Inc.
307 Westlake Ave. N., Suite 300
Seattle, Washington 98109
Telephone: 206-344-5821
Fax: 206-344-5823
www.inbios.com

Product Specification Sheet

Product	Monoclonal antibody specific to <i>Escherichia coli</i>
Catalog Number	INB9141
Description	Antibody specific to <i>Escherichia coli</i> O157 lipopolysaccharide
Specificity	LPS
Host Animal	Mouse
Format	Purified, liquid
Isotype	IgG1
Buffer	PBS pH 7.4
Application	Suitable for use in ELISA and Rapid Immunoassays. Each laboratory should determine an optimum working titer for their own application. Other formats have not been tested.
Storage	Store at -20°C. Avoid multiple freeze/thaw cycles. Prepare working dilution only prior to immediate use.

FOR RESEARCH USE ONLY. Not for use in Diagnostic Procedures.

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.